

# X-ray ICSLA\_XRY\_001

## Purpose

Construct and analyse digital X-ray images in immobilised mice using a Faxitron X-Ray system or NTB digital X-ray scanner.

## Experimental Design

- **Minimum number of animals** : 4M + 4F
- **Age at test**: Week 56

## Procedure

Perform steps 3.1 and 3.2 unless the mouse is coming from DEXA analysis

3.1 Anesthetize the mice.

3.2 Monitor the animal carefully until unconsciousness by ensuring the mouse is adequately sedated.

3.3 Mount the anesthetized mouse onto an X-ray permeable plate. Curve the tail along the left side of the mouse. (Note: this will enable identification of the X-ray as a dorsal or ventral image.)

Faxitron X-ray system

3.3.1 Insert the plate into the Faxitron X-ray system.

3.3.2 Take an image using the software that is supplied by the X-ray scanner.

3.3.3 Set the machine to a voltage of 25kV with an integration time of 40ms.

3.3.4 Take an image using the iX-Pect software that is supplied by the X-ray scanner.

NTB digital X-ray scanner

3.4 Remove the mouse once the image is captured. Place the mouse on a heated mat, set at 37°C, in a cage and monitor closely until consciousness is fully regained.

3.5 Analyze the image.

3.6 Or continue with ABR.

# Notes

The protocol described herein employs anaesthesia using intraperitoneal (IP) administration; however, anaesthesia using isoflurane is also suitable.

[Digit integrity IMPC\\_XRY\\_031\\_001](#) *Definition:* The state of being whole and undivided.

## Parameters and Metadata

### Fusion of ribs ICSLA\_XRY\_011\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: true

Options: imageOnly, not fused, unobservable, fused,

-----

### Cervical processes ICSLA\_XRY\_063\_001 | v1.1

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Options: imageOnly, normal, unobservable, abnormal,

-----

### Scanner equipment manufacturer ICSLA\_XRY\_045\_001 | v1.1

procedureMetadata

Req. Analysis: false      Req. Upload: false      Is Annotated: false

---

## Settings time of exposure ICSLA\_XRY\_039\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: s

---

## Tibia ICSLA\_XRY\_025\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: imageOnly, unobservable, normal, abnormal,

---

## Teeth ICSLA\_XRY\_005\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: imageOnly, normal, abnormal, unobservable,

---

## Number of pelvic vertebrae ICSLA\_XRY\_016\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

---

## Caudal processes ICSLA\_XRY\_067\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: abnormal, imageOnly, normal, unobservable,

---

## XRay Images Whole Body Dorso Ventral ICSLA\_XRY\_034\_001 | v1.3

seriesMediaParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Increments: Minimum 1

---

## Xray equipment model ICSLA\_XRY\_038\_001 | v1.4

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

**Options:** MX-20, MS-FX Pro, R\_mCT (Micro-CT), MX-20-DC4, UltraFocus 100, UltraFocus DXA, MX-20-DC12,

---

**Number of thoracic vertebrae** ICSLA\_XRY\_014\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** true

---

**Shape of ribs** ICSLA\_XRY\_010\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** true

**Options:** imageOnly, unobservable, normal, abnormal,

---

**Xray equipment ID** ICSLA\_XRY\_036\_001 | v1.1

procedureMetadata

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** false

---

**Number of lumbar vertebrae** ICSLA\_XRY\_015\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

---

## Comment on XRay image ICSLA\_XRY\_032\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Number of digits ICSLA\_XRY\_028\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

---

## Syndactylism ICSLA\_XRY\_029\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Options:** no syndactylism, imageOnly, Syndactylism, unobservable,

---

## Fibula ICSLA\_XRY\_026\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** normal, unobservable, imageOnly, abnormal,

---

## Clavicle ICSLA\_XRY\_007\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** unobservable, abnormal, imageOnly, normal,

---

## Thoracic vertebrae morphology ICSLA\_XRY\_071\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** abnormal, unobservable, imageOnly, normal,

---

## Lumbar vertebrae morphology ICSLA\_XRY\_072\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** abnormal, unobservable, imageOnly, normal,

---

**Alive** ICSLA\_XRY\_035\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** false

---

**Scoliosis** ICSLA\_XRY\_056\_001 | v1.2

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** true

**Options:** no, imageOnly, yes, unobservable,

---

**Free Ontology** ICSLA\_XRY\_054\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** false

---

**Lordosis** ICSLA\_XRY\_058\_001 | v1.1

simpleParameter



**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** imageOnly, unobservable, no, yes,

---

## Sacral processes ICSLA\_XRY\_066\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** normal, abnormal, unobservable, imageOnly,

---

## Scapulae ICSLA\_XRY\_006\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** normal, imageOnly, abnormal, unobservable,

---

## Pelvis ICSLA\_XRY\_012\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** imageOnly, abnormal, unobservable, normal,

---

## Polysyndactylism ICSLA\_XRY\_062\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: yes, unobservable, imageOnly, no,

---

## XRay Images Skull Dorso Ventral Orientation ICSLA\_XRY\_051\_001 | v1.2

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

---

## Tibia length ICSLA\_XRY\_033\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: mm

---

# Rib morphology ICSLA\_XRY\_069\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: unobservable, imageOnly, normal, abnormal,

---

# Radius ICSLA\_XRY\_022\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: imageOnly, unobservable, abnormal, normal,

---

# Number of cervical vertebrae ICSLA\_XRY\_013\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

---

# Caudal vertebrae morphology ICSLA\_XRY\_074\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** normal, unobservable, abnormal, imageOnly,

---

**Joints** ICSLA\_XRY\_027\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** true

**Options:** imageOnly, abnormal, unobservable, normal,

---

**Missing cranial rib** ICSLA\_XRY\_068\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** false

**Options:** Yes, No,

---

**Experimenter ID (scan)** ICSLA\_XRY\_041\_001 | v1.1

procedureMetadata

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** false

---

**XRay Images Forepaw** ICSLA\_XRY\_049\_001 | v1.1

seriesMediaParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Increments:** Minimum 1

---

## Xray equipment manufacturer ICSLA\_XRY\_037\_001 | v1.4

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** false

**Is Annotated:** false

**Options:** Faxitron X-Ray Corporation, Faxitron Bioptics LLC, Bruker, Rigaku Corporation (Micro-CT),

---

## Shape of spine ICSLA\_XRY\_055\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** imageOnly, normal, abnormal, unobservable,

---

## Fusion of vertebrae ICSLA\_XRY\_019\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** no, imageOnly, yes, unobservable,

---

**Number of ribs left** ICSLA\_XRY\_009\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

---

**Pelvic vertebrae morphology** ICSLA\_XRY\_073\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** unobservable, normal, abnormal, imageOnly,

---

**Cervical vertebrae morphology** ICSLA\_XRY\_070\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** imageOnly, unobservable, abnormal, normal,

---

# Lumbar processes ICSLA\_XRY\_065\_001 | v1.1

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Options: imageOnly, normal, abnormal, unobservable,

---

# Voltage settings ICSLA\_XRY\_043\_001 | v1.1

procedureMetadata

Req. Analysis: false      Req. Upload: true      Is Annotated: false

---

# Settings level ICSLA\_XRY\_040\_001 | v1.1

procedureMetadata

Req. Analysis: false      Req. Upload: true      Is Annotated: false

---

# Kyphosis ICSLA\_XRY\_057\_001 | v1.1

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Options: imageOnly, yes, unobservable, no,

---

**Shape of ribcage** ICSLA\_XRY\_059\_001 | v1.1

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

Options: imageOnly, normal, unobservable, abnormal,

---

**Scanner equipment ID** ICSLA\_XRY\_044\_001 | v1.1

procedureMetadata

Req. Analysis: false      Req. Upload: false      Is Annotated: false

---

**Maxilla/Pre-maxilla** ICSLA\_XRY\_003\_001 | v1.1

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: true

Options: normal, unobservable, abnormal, imageOnly,

---

**Humerus** ICSLA\_XRY\_021\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: true



**Options:** unobservable, normal, abnormal, imageOnly,

---

**Number of ribs right** ICSLA\_XRY\_008\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** true

---

**Date Xray equipment last calibrated** ICSLA\_XRY\_047\_001 | v1.3

procedureMetadata

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** false

---

**Scanner equipment model** ICSLA\_XRY\_046\_001 | v1.1

procedureMetadata

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** false

---

**XRay Images Skull Lateral Orientation** ICSLA\_XRY\_050\_001 | v1.2

seriesMediaParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** false

**Increments:** Minimum 1

---

**Brachydactyly** ICSLA\_XRY\_030\_001 | v1.1

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** imageOnly, no, yes, unobservable,

---

**Zygomatic bone** ICSLA\_XRY\_002\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** abnormal, unobservable, imageOnly, normal,

---

**Ulna** ICSLA\_XRY\_023\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** normal, unobservable, imageOnly, abnormal,

---

## Skull shape ICSLA\_XRY\_001\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: abnormal, unobservable, imageOnly, normal,

---

## Date Scanner equipment last calibrated ICSLA\_XRY\_053\_001 | v1.

1

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Thoracic processes ICSLA\_XRY\_064\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: unobservable, normal, imageOnly, abnormal,

---

## Shape of vertebrae ICSLA\_XRY\_018\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** unobservable, imageOnly, abnormal, normal,

---

## XRay Images Whole Body Lateral Orientation ICSLA\_XRY\_048\_001 | v1.2

seriesMediaParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Increments:** Minimum 1

---

## Femur ICSLA\_XRY\_024\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** unobservable, abnormal, imageOnly, normal,

---

## Digit integrity ICSLA\_XRY\_031\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** yes, unobservable, no, imageOnly,

---

## Fusion of processes ICSLA\_XRY\_061\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: yes, no, imageOnly, unobservable,

---

## Transitional vertebrae ICSLA\_XRY\_060\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: no, unobservable, imageOnly, yes,

---

## Mandibles ICSLA\_XRY\_004\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: normal, abnormal, unobservable, imageOnly,

---

# Number of caudal vertebrae ICSLA\_XRY\_017\_001 | v1.1

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: true

---

# XRay Images Hind Leg and Hip ICSLA\_XRY\_052\_001 | v1.1

seriesMediaParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

Increments: Minimum 1

---

# Processes on vertebrae ICSLA\_XRY\_020\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: true

Options: imageOnly, abnormal, unobservable, normal,

---

# Experimenter ID (analysis) ICSLA\_XRY\_042\_001 | v1.1

procedureMetadata

Req. Analysis: false      Req. Upload: true      Is Annotated: false

---

