

Eye Morphology HAS_EYE_002

Purpose

To detect abnormalities in eye morphology.

Experimental Design

- **Minimum number of animals** : 7M + 7F
- **Age at test**: Ideal age = 15 weeks \pm 3 days. Minimal age = 14 weeks
- **Sex**: We would not expect the results of this test to show sexual dimorphism

Procedure

1. Examine the anterior of both eyes (e.g. with slit lamp) and record any abnormalities
2. Test the iris/pupil light response
3. Image abnormal eyes as a minimum or all eyes if capacity permits
4. Dilate both eyes
5. Examine the anterior and posterior of both dilated eyes (e.g. with slit lamp and ophthalmoscope) and record any abnormalities
6. Image abnormal eyes as a minimum or all eyes if capacity permits

OCT:

1. Turn on the OCT and start the database
2. Anaesthetize mouse
3. Prepare mouse eyes with drops and place contact lens (focal length 10 mm) on the right eye
4. Enter mouse data in the "Create new patient file" area and switch to the "Acquisition" window
5. Move the OCT camera to the right position and activate measurement modus
6. Place mouse collaterally to the OCT camera on the right side of a platform that is fixed in front of the OCT lens
7. Search the contact lens in the live picture of the fundus image field and place the pupil of the mouse eye in the centre of the window
8. Move the OCT camera such that OCT lens and contact lens touch each other
9. Focus the fundus picture by slightly moving up/down or forward/backward
10. Save fundus images
11. Set the „Ref.Arm“ ruler such that the section of the retina is placed in the centre of the blue rectangle
12. Set the mode of measurement on „vertical, horizontal line“
13. Move the blue horizontal line in the fundus image field to the optic nerve level
14. Save images of retinal sections
15. Move the OCT camera to the left position
16. Repeat measurement procedure for the left eye

Scheimpflug Imaging:

1. Turn on the Pentacam and start the patient data management
2. Apply one drop 0.5% Atropine to each mouse eye for pupil dilation
3. Enter mouse data in the "Patient" group box and switch to the Scan menu
4. Activate the "1 Picture" modus in the "Image Options" area
5. Move Pentacam to the right position
6. Hold the mouse on a platform such that the vertical LED 475 nm light slit is orientated in the center of the right eye ball
7. Guarantee optimal focus by using the fine adjustment software tool in the adjustment window
8. Start imaging manually by pressing the "Start Scan" button
9. Scheimpflug images are saved automatically
10. Move Pentacam to the left position
11. Repeat measurement procedure for the left eye

Notes

- As a minimum, all abnormalities should be imaged.
 - Where capacity permits, all mice can be imaged
- Majority of parameters can be analysed using the standard approach for assessing categorical data. To increase power for analysis purposes, where an abnormality is detected in the left, right or both eyes, the data may be combined to generate one "abnormal" category.

Data QC

Image QC is typically performed during data collection to ensure high quality images are captured whilst eyes are dilated etc.

Parameters and Metadata

Eye HAS_EYE_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: eye

Options: present, absent left eye, absent right eye, absent both eyes,

Bulging eye HAS_EYE_002_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: bulging_eye

Options: absent, no data left eye, no data right eye, present left eye, present right eye, present both eyes, no data for both eyes,

Eye Hemorrhage or Blood Presence HAS_EYE_003_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: eye_hemorrhage_or_blood_presence

Options: absent, no data left eye, no data right eye, present left eye, present right eye, present both eyes, no data for both eyes,

Eyelid morphology HAS_EYE_004_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: eyelid_morphology

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Eyelid closure HAS_EYE_005_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: eyelid_closure

Options: normal, no data left eye, no data right eye, left eye closed, right eye closed, both eyes closed, no data for both eyes,

Narrow eye opening HAS_EYE_006_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: narrow_eye_opening

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Cornea HAS_EYE_007_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: cornea

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Corneal opacity HAS_EYE_008_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: corneal_opacity

Options: absent, no data left eye, no data right eye, present left eye, present right eye, present both eyes, no data for both eyes,

Corneal vascularization HAS_EYE_009_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: corneal_vascularization

Options: absent, no data left eye, no data right eye, present left eye, present right eye, present both eyes, no data for both eyes,

Iris/Pupil HAS_EYE_010_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: iris_pupil

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Pupil Position HAS_EYE_011_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: pupil_position

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Pupil Shape HAS_EYE_012_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: pupil_shape

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Pupil Dilation HAS_EYE_013_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: pupil_dilation

Options: normal, no data left eye, no data right eye, left eye dilated, right eye dilated, both eyes dilated, no data for both eyes,

Pupil Light Response HAS_EYE_014_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: pupil_light_response

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Iris Pigmentation HAS_EYE_015_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: iris_pigmentation

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Lens HAS_EYE_016_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: lens

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Lens Opacity HAS_EYE_017_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: lens_opacity

Options: absent, no data left eye, no data right eye, present left eye, present right eye, present both eyes, no data for both eyes,

Fusion between cornea and lens HAS_EYE_018_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: fusion_between_cornea_and_lens

Options: absent, no data left eye, no data right eye, present left eye, present right eye, present both eyes, no data for both eyes,

Synechia HAS_EYE_019_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: synechia

Options: absent, no data left eye, no data right eye, present left eye, present right eye, present both eyes, no data for both eyes,

Retina HAS_EYE_020_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: retina

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Retinal Pigmentation HAS_EYE_021_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: retinal_pigmentation

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Retinal Structure HAS_EYE_022_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: retinal_structure

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Optic Disc HAS_EYE_023_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: optic_disc

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Retinal Blood Vessels HAS_EYE_024_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: retinal_blood_vessels

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Retinal Blood Vessels Structure HAS_EYE_025_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: retinal_blood_vessels_structure

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Retinal Blood Vessels Pattern HAS_EYE_026_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: retinal_blood_vessels_pattern

Options: normal, no data left eye, no data right eye, left eye abnormal, right eye abnormal, both eyes abnormal, no data for both eyes,

Persistence of hyaloid vascular system HAS_EYE_027_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: persistence_of_hyaloid_vascular_system

Options: absent, no data left eye, no data right eye, present left eye, present right eye, present both eyes, no data for both eyes,

Slit Lamp observation HAS_EYE_028_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: slit_lamp_observation

Ophthalmoscope Observation HAS_EYE_029_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ophthalmoscope_observation

Slit Lamp Equipment ID HAS_EYE_030_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: slit_lamp_equipment_id

Slit Lamp Equipment Manufacturer HAS_EYE_031_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: slit_lamp_equipment_manufacturer

Slit Lamp Equipment Model HAS_EYE_032_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: slit_lamp_equipment_model

Ophthalmoscope Equipment ID HAS_EYE_033_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ophthalmoscope_equipment_id

Ophthalmoscope Equipment Manufacturer HAS_EYE_034_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ophthalmoscope_equipment_manufacturer

Ophthalmoscope Equipment Model HAS_EYE_035_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ophthalmoscope_equipment_model

PIL number HAS_EYE_036_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: experimenter_id

Optical Coherence Tomography Equipment ID HAS_EYE_037_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: optical_coherence_tomography_equipment_id

Optical Coherence Tomography Equipment Manufacturer HAS_EYE_038_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: optical_coherence_tomography_equipment_manufacturer

Optical Coherence Tomography Equipment Model HAS_EYE_039_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: optical_coherence_tomography_equipment_model

Scheimpflug Equipment ID HAS_EYE_040_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: scheimpflug_equipment_id

Scheimpflug Equipment Manufacturer HAS_EYE_041_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: scheimpflug_equipment_manufacturer

Scheimpflug Equipment Model HAS_EYE_042_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: scheimpflug_equipment_model

Dilation Method HAS_EYE_043_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: dilation_method

Topical Anesthetic HAS_EYE_044_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: topical_anesthetic

General Anesthetic HAS_EYE_045_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: general_anesthetic

Date Ophthalmoscope equipment last calibrated HAS_EYE_047_001 | v1.0
procedureMetadata

Req. Analysis: false Req. Upload: false Is Annotated: false

Date Scheimpflug equipment last calibrated HAS_EYE_048_001 | v1.0
procedureMetadata

Req. Analysis: false Req. Upload: false Is Annotated: false

Date OCT equipment last calibrated HAS_EYE_049_001 | v1.0
procedureMetadata

Req. Analysis: false Req. Upload: false Is Annotated: false

Images Ophthalmoscopy HAS_EYE_050_001 | v1.0
seriesMediaParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Increments: Minimum 1

Images Slit Lamp HAS_EYE_051_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

Sheimpflug Lens description HAS_EYE_052_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Scheimpflug description HAS_EYE_053_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Scheimpflug min left eye lens density HAS_EYE_054_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

Scheimpflug max left eye lens density HAS_EYE_055_001 | v1.0
simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: %

Scheimpflug mean left eye lens density HAS_EYE_056_001 | v1.0
simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: %

Scheimpflug min right eye lens density HAS_EYE_057_001 | v1.0
simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: %

Scheimpflug max right eye lens density HAS_EYE_058_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

Scheimpflug mean right eye lens density HAS_EYE_059_001 | v1.

0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

OCT right corneal thickness HAS_EYE_060_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT right anterior chamber depth HAS_EYE_061_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT right total retinal thickness HAS_EYE_062_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT right inner nuclear layer HAS_EYE_063_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT right outer nuclear layer HAS_EYE_064_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT right posterior chamber depth HAS_EYE_065_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT left corneal thickness HAS_EYE_066_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT left anterior chamber depth HAS_EYE_067_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT left total retinal thickness HAS_EYE_068_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT left inner nuclear layer HAS_EYE_069_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT left outer nuclear layer HAS_EYE_070_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT left posterior chamber depth HAS_EYE_071_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

OCT B-scan of right retina HAS_EYE_072_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

OCT B-scan of left retina HAS_EYE_073_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

OCT VIP of right fundus HAS_EYE_074_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

OCT VIP of left fundus HAS_EYE_075_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

OCT B-scan of right cornea and lens HAS_EYE_076_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

OCT B-scan of left cornea and lens HAS_EYE_077_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

OCT VIP of right eye HAS_EYE_078_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

OCT VIP of left eye HAS_EYE_079_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

Corneal Sclerization HAS_EYE_080_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: absent, no data left eye, no data right eye, no data for both eyes, present left eye, present right eye, present both eyes,

Corneal deposits HAS_EYE_081_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: absent, no data left eye, no data right eye, no data for both eyes, present left eye, present right eye, present both eyes,

Iris transillumination HAS_EYE_082_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: normal, no data left eye, no data right eye, no data for both eyes, left eye abnormal, right eye abnormal, both eyes abnormal,

Vitreous HAS_EYE_083_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: normal, no data left eye, no data right eye, no data for both eyes, left eye abnormal, right eye abnormal, both eyes abnormal,

Date of procedure HAS_EYE_046_001 | v1.1

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

General comments HAS_EYE_084_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Procedural comments HAS_EYE_085_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: false