

# Body Weight HAS\_BWT\_002

## Purpose

The body weight test measures the weight of the mouse in a time series, allowing monitoring of its evolution; also, it is required in many other procedures.

## Equipment

- Laboratory Balance

## Procedure

- Place a cage containing mice on an operation table
- Remove a mouse from its cage. Determine and record its weight.
- Return the mouse to the cage or continue with a procedure.

## Notes

The body weight of each mouse is to be measured weekly, from the 4<sup>th</sup> until the 16<sup>th</sup> week. Optional additional weighs may also be recorded.

## Parameters and Metadata

### Body weight HAS\_BWT\_001\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: g

Description: body\_weight

---

### General comments about the mouse HAS\_BWT\_002\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Description:** general\_comments\_about\_the\_mouse

---

## Equipment ID HAS\_BWT\_003\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Description:** equipment\_name

---

## Equipment manufacturer HAS\_BWT\_004\_001 | v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** equipment\_manufacturer

**Options:** Ohaus, Denver Instrument, Scienteck, Mettler Toledo,

---

## Experimenter ID HAS\_BWT\_005\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: experimenter\_id

Date equipment last calibrated

HAS\_BWT\_006\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Equipment model

HAS\_BWT\_007\_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Options: AV2101, SI-2002, SL-3100D, TP-202, PG3001-S,

Date of procedure

HAS\_BWT\_008\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Procedural comments

HAS\_BWT\_009\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

-----