

Eye Morphology RBRCLA_EYE_001

Purpose

To detect abnormalities in eye morphology.

Experimental Design

- **Minimum number of animals :** 4M + 4F
- **Age at test:** Week 53
- **Sex:** We do not expect the results of this test to show sexual dimorphism

Procedure

1. Examine the anterior of both eyes (e.g. with slit lamp) and record any abnormalities
2. Test the iris/pupil light response
3. Image abnormal eyes as a minimum or all eyes if capacity permits
4. Dilate both eyes
5. Examine the anterior and posterior of both dilated eyes (e.g. with slit lamp and ophthalmoscope) and record any abnormalities
6. Image abnormal eyes as a minimum or all eyes if capacity permits

OCT:

1. Turn on the OCT and start the database
2. Anaesthetize mouse
3. Prepare mouse eyes with drops and place contact lens (focal length 10 mm) on the right eye
4. Enter mouse data in the "Create new patient file" area and switch to the "Acquisition" window
5. Move the OCT camera to the right position and activate measurement mode
6. Place mouse collaterally to the OCT camera on the right side of a platform that is fixed in front of the OCT lens
7. Search the contact lens in the live picture of the fundus image field and place the pupil of the mouse eye in the centre of the window
8. Move the OCT camera such that OCT lens and contact lens touch each other
9. Focus the fundus picture by slightly moving up/down or forward/backward
10. Save fundus images
11. Set the "Ref.Arm" ruler such that the section of the retina is placed in the centre of the blue rectangle
12. Set the mode of measurement on "vertical, horizontal line"
13. Move the blue horizontal line in the fundus image field to the optic nerve level
14. Save images of retinal sections
15. Move the OCT camera to the left position

16. Repeat measurement procedure for the left eye

Scheimpflug Imaging:

1. Turn on the Pentacam and start the patient data management
2. Apply one drop 0.5% Atropine to each mouse eye for pupil dilation
3. Enter mouse data in the "Patient" group box and switch to the Scan menu
4. Activate the "1 Picture" modus in the "Image Options" area
5. Move Pentacam to the right position
6. Hold the mouse on a platform such that the vertical LED 475 nm light slit is orientated in the center of the right eye ball
7. Guarantee optimal focus by using the fine adjustment software tool in the adjustment window
8. Start imaging manually by pressing the "Start Scan" button
9. Scheimpflug images are saved automatically
10. Move Pentacam to the left position
11. Repeat measurement procedure for the left eye

Notes

- Only the eye that was not used for retro-orbital bleeds at early adult should be examined.
- As a minimum, all abnormalities should be imaged.
 - Where capacity permits, all mice can be imaged
- Majority of parameters can be analysed using the standard approach for assessing categorical data. To increase power for analysis purposes, where an abnormality is detected in the left, right or both eyes, the data may be combined to generate one "abnormal" category.

Data QC

Image QC is typically performed during data collection to ensure high quality images are captured whilst eyes are dilated etc.

Parameters and Metadata

Eye RBRCLA_EYE_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: eye

Options: present, absent left eye, absent right eye, absent both eyes,

Bulging eye RBRCLA_EYE_002_001 | v1.0
simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: bulging_eye

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Eye Hemorrhage or Blood Presence RBRCLA_EYE_003_001 | v1.0
simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: eye_hemorrhage_or_blood_presence

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Eyelid morphology RBRCLA_EYE_004_001 | v1.0
simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: eyelid_morphology

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Eyelid closure RBRCLA_EYE_005_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: eyelid_closure

Options: normal, left eye closed, right eye closed, both eyes closed, no data left eye,
no data right eye, no data for both eyes, no data left eye, right eye closed,
no data right eye, left eye closed,

Narrow eye opening RBRCLA_EYE_006_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: narrow_eye_opening

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Cornea RBRCLA_EYE_007_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: cornea

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Corneal opacity RBRCLA_EYE_008_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: corneal_opacity

Options: absent, present left eye, present right eye, present both eyes, no data left eye,
no data right eye, no data for both eyes, no data left eye, present right eye,
no data right eye, present left eye,

Corneal vascularization RBRCLA_EYE_009_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: corneal_vascularization

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Iris/Pupil RBRCLA_EYE_010_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: iris_pupil

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal, no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal, no data right eye, left eye abnormal,

Pupil Position RBRCLA_EYE_011_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: pupil_position

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal, no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal, no data right eye, left eye abnormal,

Pupil Shape RBRCLA_EYE_012_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: pupil_shape

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Pupil Dilation RBRCLA_EYE_013_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: pupil_dilation

Options: normal, left eye dilated, right eye dilated, both eyes dilated, no data left eye,
no data right eye, no data for both eyes, no data left eye, right eye dilated,
no data right eye, left eye dilated,

Pupil Light Response RBRCLA_EYE_014_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: pupil_light_response

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,

no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,

no data right eye, left eye abnormal,

Iris Pigmentation RBRCLA_EYE_015_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: iris_pigmentation

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,

no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,

no data right eye, left eye abnormal,

Lens RBRCLA_EYE_016_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: lens

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,

no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,

no data right eye, left eye abnormal,

Lens Opacity RBRCLA_EYE_017_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: lens_opacity

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Fusion between cornea and lens RBRCLA_EYE_018_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: fusion_between_cornea_and_lens

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Synechia RBRCLA_EYE_019_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: synechia

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Retina RBRCLA_EYE_020_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: retina

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal, no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal, no data right eye, left eye abnormal,

Retinal Pigmentation RBRCLA_EYE_021_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: retinal_pigmentation

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal, no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal, no data right eye, left eye abnormal,

Retinal Structure RBRCLA_EYE_022_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Description: retinal_structure

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Optic Disc RBRCLA_EYE_023_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Description: optic_disc

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Retinal Blood Vessels RBRCLA_EYE_024_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Description: retinal_blood_vessels

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Retinal Blood Vessels Structure RBRCLA_EYE_025_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Description: retinal_blood_vessels_structure

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Retinal Blood Vessels Pattern RBRCLA_EYE_026_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: retinal_blood_vessels_pattern

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Persistence of hyaloid vascular system RBRCLA_EYE_027_001 |

v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: persistence_of_hyaloid_vascular_system

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Slit Lamp observation RBRCLA_EYE_028_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: slit_lamp_observation

Ophthalmoscope Observation RBRCLA_EYE_029_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ophthalmoscope_observation

Slit Lamp Equipment ID RBRCLA_EYE_030_001 | v1.0

procedureMetadata

Req. Analysis: false Req. Upload: false Is Annotated: false

Description: slit_lamp_equipment_id

Slit Lamp Equipment Manufacturer RBRCLA_EYE_031_001 | v1.0

procedureMetadata

Req. Analysis: true Req. Upload: false Is Annotated: false

Description: slit_lamp_equipment_manufacturer

Options: Zeiss, Haag-Streit, MuLe, Kowa, CSO, Phoenix Research Labs, Topcon,

Slit Lamp Equipment Model RBRCLA_EYE_032_001 | v1.0

procedureMetadata

Req. Analysis: true Req. Upload: false Is Annotated: false

Description: slit_lamp_equipment_model

Options: SL30, SL130, BQ 900 LED/IM-900, S350, SL-15, SL 990, SL 139, 30 SL-M, Micron III slit lamp extension, SL-7E,

Ophthalmoscope Equipment ID RBRCLA_EYE_033_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ophthalmoscope_equipment_id

Ophthalmoscope Equipment Manufacturer RBRCLA_EYE_034_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Description: ophthalmoscope_equipment_manufacturer

Options: Haag-Streit, Heine, Phoenix, Kowa, Karl Storz / Nikon, Phoenix Research Labs, Heine / Volk, Keeler LTD,

Ophthalmoscope Equipment Model RBRCLA_EYE_035_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Description: ophthalmoscope_equipment_model

Options: Sigma 150K, Omega 500 Unplugged, Micron III, Genesis-D, OMEGA 180 / Superfield NC, Xenon Nova 175W light source + HOPKINS optic 1218AA /Nikon D5100 + 85 mm f/1.8 lens, Omega 180 / 60D, SL4 4AA, Genesis, Genesis-DF, Micron IV,

Experimenter ID RBRCLA_EYE_036_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: experimenter_id

Optical Coherence Tomography Equipment ID RBRCLA_EYE_037_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: optical_coherence_tomography_equipment_id

Optical Coherence Tomography Equipment Manufacturer RBRCLA_EYE_038_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Description: optical_coherence_tomography_equipment_manufacturer

Options: Bioptigen, Heidelberg Engineering,

Optical Coherence Tomography Equipment Model

RBRCLA_EYE_039_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Description: optical_coherence_tomography_equipment_model

Options: EnvisuTM R-Series SDOIS, Envisu R2200, Spectralis,

Scheimpflug Equipment ID

RBRCLA_EYE_040_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: scheimpflug_equipment_id

Scheimpflug Equipment Manufacturer

RBRCLA_EYE_041_001 | v1

.0

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Description: scheimpflug_equipment_manufacturer

Options: Oculus GmbH,

Scheimpflug Equipment Model RBRCLA_EYE_042_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Description: scheimpflug_equipment_model

Options: Pentacam,

Dilation Method RBRCLA_EYE_043_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: dilation_method

Options: Atropine, Tropicamide, Tropicamide+Phenylephrin, Cyclopentolate hydrochloride, Phenylephrine hydrochloride, None, Atropine sulphate, Cyclopentolate hydrochloride+Phenylephrine hydrochloride,

Topical Anesthetic RBRCLA_EYE_044_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: topical_anesthetic

Options: Atropine, Oxybuprocain, Mydriacyl, Phenylephrine hydrochloride, Hydrochloride, No anesthesia, Atropine sulphate,

General Anesthetic RBRCLA_EYE_045_001 | v1.0

procedureMetadata

Req. Analysis: true **Req. Upload:** true **Is Annotated:** false

Description: general_anesthetic

Options: Ketamine+Xylazine, Isoflurane, Euthatal, Avertin, No anesthesia, Ketamine+Medetomidine,

Date Slit Lamp equipment last calibrated RBRCLA_EYE_046_001

| v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Date Ophthalmoscope equipment last calibrated RBRCLA_EY

E_047_001 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Date Scheimpflug equipment last calibrated RBRCLA_EYE_048_

001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Date OCT equipment last calibrated RBRCLA_EYE_049_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Images Ophthalmoscopy RBRCLA_EYE_050_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

Images Slit Lamp RBRCLA_EYE_051_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

Sheimpflug Lens description

RBRCLA_EYE_052_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Scheimpflug description

RBRCLA_EYE_053_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Min left eye lens density

RBRCLA_EYE_054_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

Max left eye lens density

RBRCLA_EYE_055_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

Mean left eye lens density RBRCLA_EYE_056_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

Min right eye lens density RBRCLA_EYE_057_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

Max right eye lens density RBRCLA_EYE_058_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

Mean right eye lens density RBRCLA_EYE_059_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

Right corneal thickness RBRCLA_EYE_060_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

Right anterior chamber depth RBRCLA_EYE_061_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

Right total retinal thickness RBRCLA_EYE_062_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: um

Right inner nuclear layer RBRCLA_EYE_063_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: um

Right outer nuclear layer RBRCLA_EYE_064_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: um

Right posterior chamber depth RBRCLA_EYE_065_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

Left corneal thickness RBRCLA_EYE_066_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

Left anterior chamber depth RBRCLA_EYE_067_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

Left total retinal thickness RBRCLA_EYE_068_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

Left inner nuclear layer RBRCLA_EYE_069_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

Left outer nuclear layer RBRCLA_EYE_070_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

Left posterior chamber depth RBRCLA_EYE_071_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

B-scan of right retina RBRCLA_EYE_072_001 | v1.0

seriesMediaParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Increments: Minimum 1

B-scan of left retina RBRCLA_EYE_073_001 | v1.0

seriesMediaParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Increments: Minimum 1

VIP of right fundus RBRCLA_EYE_074_001 | v1.0

seriesMediaParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Increments: Minimum 1

VIP of left fundus RBRCLA_EYE_075_001 | v1.0

seriesMediaParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Increments: Minimum 1

B-scan of right cornea and lens RBRCLA_EYE_076_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

B-scan of left cornea and lens RBRCLA_EYE_077_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

VIP of right eye RBRCLA_EYE_078_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

VIP of left eye RBRCLA_EYE_079_001 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Increments: Minimum 1

Corneal Sclerization RBRCLA_EYE_080_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Corneal deposits RBRCLA_EYE_081_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Iris transillumination RBRCLA_EYE_082_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Vitreous RBRCLA_EYE_083_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: normal, left eye abnormal, right eye abnormal, both eyes abnormal,
no data left eye, no data right eye, no data for both eyes, no data left eye, right eye abnormal,
no data right eye, left eye abnormal,

Corneal mineralization RBRCLA_EYE_084_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: absent, present left eye, present right eye, present both eyes, no data left eye,
no data right eye, no data for both eyes, no data left eye, present right eye,
no data right eye, present left eye,

Corneal ulcer RBRCLA_EYE_085_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Lacrimation RBRCLA_EYE_086_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: absent, present left eye, present right eye, present both eyes, no data left eye, no data right eye, no data for both eyes, no data left eye, present right eye, no data right eye, present left eye,

Right vitreous humor thickness RBRCLA_EYE_087_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: um

Left vitreous humour thickness RBRCLA_EYE_088_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: um

Ophthalmoscope Lens Model RBRCLA_EYE_089_001 | v1.0

procedureMetadata

Req. Analysis: false Req. Upload: false Is Annotated: false

Right eye diameter RBRCLA_EYE_090_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: mm

Left eye diameter RBRCLA_EYE_091_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: mm

Retina (combined) RBRCLA_EYE_092_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Derivation:

retinaCombined('RBRCLA_EYE_020_001', 'RBRCLA_EYE_021_001',
'RBRCLA_EYE_022_001')
