

Ophthalmoscope M-G-P_013_001

Parameters and Metadata

Eye M-G-P_013_001_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: true

Description: Eye

Options: normal, abnormal,

Description M-G-P_013_001_002 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Description: Description

Left Retina M-G-P_013_001_003 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Description: Left_Retina

Options: normal, abnormal, no data,

Left Retina Pigmentation M-G-P_013_001_004 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Retina_Pigmentation

Options: normal, abnormal, not defined, no data,

Left Retina Structure M-G-P_013_001_005 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Retina_Structure

Options: normal, degenerate, not defined, no data,

Left Optic Disc M-G-P_013_001_006 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Optic_Disc

Options: normal, abnormal, not defined, no data,

Right Retina M-G-P_013_001_007 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Retina

Options: normal, abnormal, no data,

Right Retina Pigmentation M-G-P_013_001_008 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Retina_Pigmentation

Options: normal, abnormal, not defined, no data,

Right Retina Structure M-G-P_013_001_009 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Retina_Structure

Options: normal, degenerate, not defined, no data,

Right Optic Disc M-G-P_013_001_010 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Optic_Disc

Options: normal, abnormal, not defined, no data,

Left Eye Blood Vessels M-G-P_013_001_011 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Eye_Blood_Vessels

Options: normal, abnormal, no data,

Left Eye Blood Vessels Structure M-G-P_013_001_012 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Eye_Blood_Vessels_Structure

Options: normal appearance, abnormal appearance, not defined, no data,

Left Eye Blood Vessels Pattern M-G-P_013_001_013 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Left_Eye_Blood_Vessels_Pattern

Options: normal, abnormal, not defined, no data,

Right Eye Blood Vessels M-G-P_013_001_014 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Eye_Blood_Vessels

Options: normal, abnormal, no data,

Right Eye Blood Vessels Structure M-G-P_013_001_015 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Eye_Blood_Vessels_Structure

Options: normal appearance, abnormal appearance, not defined, no data,

Right Eye Blood Vessels Pattern M-G-P_013_001_016 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Right_Eye_Blood_Vessels_Pattern

Options: normal, abnormal, not defined, no data,

Digital Image/s M-G-P_013_001_017 | v1.0

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: Digital_Image/s

Body Weight M-G-P_013_001_018 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: g

Description: Body_Weight

Equipment name M-G-P_013_001_801 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_name

Equipment manufacturer M-G-P_013_001_802 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_manufacturer

Equipment model M-G-P_013_001_803 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_model

Topical agent(s) M-G-P_013_001_804 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: Topical_agent_s_

Options: atropine, tropicamide, atropine + neosynephrine, tropicamide + neosynephrine, tropicamide + phenylephrine,
