

DEXA HAROLD_005_001

Parameters and Metadata

Body weight HAROLD_005_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

Description: Weight

Fat mass HAROLD_005_001_002 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: g

Description: Fat_Mass

Lean mass HAROLD_005_001_003 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: g

Description: Lean_Mass

Bone Mineral Density (excluding skull) HAROLD_005_001_004 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: g/cm^2

Description: Bone_Mineral_Density_Excluding_skull

Bone Mineral Content HAROLD_005_001_005 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: g

Description: Bone_Mineral_Content

Body length HAROLD_005_001_006 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: cm

Description: Body_Length

Date of procedure HAROLD_005_001_007 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: date_of_procedure

BMC/Body weight HAROLD_005_001_701 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: BMC/Weight

Derivation: archived('(HAROLD_005_001_005/HAROLD_005_001_001)')

Lean/Body weight HAROLD_005_001_702 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Lean/Weight

Derivation: archived('(HAROLD_005_001_003/HAROLD_005_001_001)')

Fat/Body weight HAROLD_005_001_703 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: Fat/Weigth

Derivation: archived('(HAROLD_005_001_002/HAROLD_005_001_001)')

Bone area (BMC/BMD) HAROLD_005_001_704 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Description: area=bmc/bmd

Derivation: archived('(HAROLD_005_001_005/HAROLD_005_001_004)')

Equipment name HAROLD_005_001_801 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_name

Equipment manufacturer HAROLD_005_001_802 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_manufacturer

Equipment model HAROLD_005_001_803 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Equipment_model

Mouse status HAROLD_005_001_804 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: Mouse_Status

Options: anaesthetised, dead,

.....