

# Hot Plate ICSLA\_HOT\_001

## Parameters and Metadata

### Type of recording ICSLA\_HOT\_008\_001 | v1.0

procedureMetadata

Req. Analysis: true                      Req. Upload: false                      Is Annotated: false

Options: Manual,

---

### Hot plate latency ICSLA\_HOT\_003\_001 | v1.0

simpleParameter

Req. Analysis: false                      Req. Upload: true                      Is Annotated: false

Unit Measured: s

---

### Type of response ICSLA\_HOT\_004\_001 | v1.0

simpleParameter

Req. Analysis: false                      Req. Upload: true                      Is Annotated: false

Options: jump, shake, lick, other,

---

## Equipment model ICSLA\_HOT\_006\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Options: LE 7406,

---

## Equipment manufacturer ICSLA\_HOT\_005\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Options: Bioseb,

---

## Experimenter ID ICSLA\_HOT\_001\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Body weight ICSLA\_HOT\_002\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

-----

**Plate temperature** ICSLA\_HOT\_009\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: C

Options: 52,

-----