

Sleep: light pulse HAS_PLS_001

Parameters and Metadata

Date of procedure HAS_PLS_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Coffin number HAS_PLS_002_001 | v1.1

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Options: 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20,

Apparatus HAS_PLS_003_001 | v1.1

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Distance travelled HAS_PLS_004_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: cm

Description: distance_travelled

Mean speed HAS_PLS_005_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Time mobile HAS_PLS_006_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: s

Time immobile HAS_PLS_007_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: s

Mobile episodes HAS_PLS_008_001 | v1.3

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Immobile episodes HAS_PLS_009_001 | v1.3

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Time segment of test HAS_PLS_010_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: s

Options: 0 - 600, 600 - 1200, 1200 - 1800, 1800 - 2400, 2400 - 3000, 3000 - 3600, 3600 - 4200, 4200 - 4800, 4800 - 5400, 5400 - 6000, 6000 - 6600, 6600 - 7200, 7200 - 7800, 7800 - 8400, 8400 - 9000, 9000 - 9600, 9600 - 10200, 10200 - 10800, 10800 - 11400, 11400 - 12000, 12000 - 12600,

General comments about the mouse HAS_PLS_011_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: general_comments_about_the_mouse

Equipment manufacturer HAS_PLS_012_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: equipment_manufacturer

Equipment ID HAS_PLS_013_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: equipment_name

Equipment model HAS_PLS_014_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Procedural comments HAS_PLS_015_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

? HAS_PLS_016_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: Pre DP, Post DP, LP,
