

# Hot Plate ICS\_HOT\_001

## Parameters and Metadata

### Body weight ICS\_HOT\_002\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: false

Unit Measured: g

---

### Hot plate latency ICS\_HOT\_003\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: false

Unit Measured: s

---

### Type of response ICS\_HOT\_004\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: false

Options: shake, lick, jump, other,

---

## Equipment manufacturer ICS\_HOT\_005\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Options: Bioseb,

---

## Equipment model ICS\_HOT\_006\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Options: LE 7406,

---

## Type of recording ICS\_HOT\_008\_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Options: Manual,

---

## Plate temperature ICS\_HOT\_009\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: C

Options: 52,

-----

**Experimenter ID** ICS\_HOT\_001\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

-----